


X-RAY Touchscreen

Control Console



A touch of the screen is all it takes.



The newly-designed CPI Touchscreen comes with all these great features as standard.

- Settings and controls right at your fingertips
- No moving parts for quiet reliable operation
- Flash memory - no hard drive
- Same great user interface with more power and speed
- Embedded GenWare™ Service Software
- Total console control and viewing angle flexibility
- Larger/brighter screen display (15") in a light compact design

All in one package!

100® & CMP200® family of X-Ray 'Generators-of-Choice'.

**Built-in functionality optimizes set up of your Touchscreen.
Optional features allow for customization.**



Anatomical Program (APR) Editor

Allows you to select, edit and store all your anatomical settings.

Customization

Create your organization's 'look-and-feel' by adding your own company name, logo, and buttons (shape & colour).



Integration

Another customized option that allows you to program your Touchscreen to provide system integration and control.



APR Language Selection

English, Italian, French, Spanish, Swedish, Chinese, Russian and more ... all at your fingertips.



On-Screen Multi-Language Keyboard



System Utilities Screen



Auto Positioning Set up Screen

**X-RAY
Generators
of Choice**

- **Easy to use**
- **Flexible**
- **Powerful**



Handswitch
(Optional)




Wall Mount (Optional)



The Indico 100® power supply
‘The heart of the CPI system’

Interface Connectivity



All dimensions are in mm (inches)



Communications & Power Industries Canada Inc.
communications & medical products division
45 River Drive / Georgetown, Ontario / Canada L7G 2J4
Hot Line Telephone: 1-888-CPI-XRAY
TEL: 905-877-0161 / FAX: 905-877-5327
E-MAIL: marketing@cmp.cpii.com
WEB: www.cpii.com/cmp

